

Department of Laboratory Services Forms

Contents	Pg
Herpesvirus Antigen Detection by Fluorescence	2
Pertussis (Whooping Cough) Agent Detection by Fluorescence	2
Group A Streptococcus Agent Detection by Fluorescence	3
Viral & Rickettsial Serology To Detect Antibody	3
Viral Isolation	4
Rabies Agent Detection And/Or Isolation To Identify	
Infection In Animals	7
Human Immunodeficiency Virus (HIV) Serology	9
Prenatal Profile	11
HBsAg (Hepatitis B surface antigen), Anti-HBs (Antibody to HbsAg), or Anti-HBc (Antibody to HB core antigen)	13
ABO , Rh factor, Antibody Screen	13
Syphilis Testing	14
Dental Fluoride (Supplement Program)	17
Smear exam for Neisseria gonorrhoeae (GC)	19
Salmonella/Shigella species and other enteric pathogens	19
Intestinal Parasites Examination (See Notes)	20
Pinworms	21
Mycobacteriology (TB) Smear & Culture	24
Chlamydia trachomatis & Neisseria gonorrhoeae	26
Newborn Screening	29
Sickle cell Screening	30
Cholesterol Screening (Must specify Total Cholesterol or Lipid Profile)	32
Lithium Level	33
Glucose Tolerance , Prenatal / Postpartum, One hour post 50 gram load Glucose	34
Fasting / Random Glucose	35

Herpesvirus Antigen Detection by Fluorescence

Patients Qualifying:

- All Local Health Department patients

Specimen:

- Slide prepared from vesicular fluid collection

Collection Kit (Herpes FA Slide Kit) Furnished by State Lab Contains:

- Microscope slide
- Styrofoam slide box
- Lab Form # 275 Viral Isolation & Immunology
- Address label # 356
- Instruction Sheet

Pertussis (Whooping Cough) Agent Detection by Fluorescence

Patients Qualifying:

- All Local Health Department patients

Specimen:

- Slide prepared from nasopharyngeal swab specimen

Collection Kit (Pertusis FA Slide Kit) Furnished by State Lab Contains:

- Microscope slide
- Styrofoam slide box
- Lab Form # 275 Viral Isolation & Immunology
- Address label # 356
- Instruction Sheet

Group A Streptococcus Agent Detection by Fluorescence

Patients Qualifying:

- All Local Health Department patients

Specimen:

- Throat Swab (transported in silica gel)

Collection Kit (Group A Strep Kit) Furnished by State Lab Contains:

- Sterile Swab for specimen collection with silica gel preservative packet
- 9" X 7" Manila Envelope
- Lab Form # 275- Viral Isolation & Immunology
- Address label # 356
- Instruction Sheet

Viral & Rickettsial Serology To Detect Antibody (see alphabetical "Reference List" for IgG & IgM Serology testing available)

Patients Qualifying:

- All Local Health Department patients

Specimen:

- Serum or Whole Blood

Collection Kit (Viral & Rickettsial serology Kit) Furnished by State Lab Contains:

- Sterile Red-stopper tube
- Lab Form # 275- Viral Isolation & Immunology
- Address label # 356
- 2-part Mailing Canister -inner & outer canisters
- Absorbent material to be provided by submitter

Viral Isolation

Patients Qualifying:

- All Local Health Department patients

Specimen:

- Throat Swab, Rectal Swab, Spinal Fluid or Tissue as appropriate to clinical symptoms

Collection Kit (Virus isolation swab) Furnished by State Lab Contains:

- Sterile Viral Isolation Swab collection outfit *
 - Lab Form # 275- Viral Isolation & Immunology (Instructions on back)
 - Address label # 356
 - Styrofoam Mailer 8" x 5" x 6.5"
 - Cold Pack
- * for non-swab specimens submitter must provide a sterile leak-proof container

Lab 275 (Rev 4/00)		Kentucky Public Health Laboratory	
VIRAL ISOLATION & IMMUNOLOGY		North Loading Dock P O Box 2020 Frankfort, KY 40601-2020 502/564-4446 Samuel B. Gregorio, Dr. P.H., Director For hand delivery : 100 Sower Blvd Suite 204	
(Submit white copy with specimen)			
Patient Information: (can use label with complete info)			
Name (Last, First, MI)		please print	
Social Security #	Sex	EO	Age (dd-mm-yyy)
Home Address			
City			
State	ZIP	County	
Send Reports to: Submitter Name or Facility			
Street Address / P O Box			
City			
State	ZIP		
Attending Physician's Name (if other than Submitter)			
EXAMINATION REQUESTED		SPECIMEN REQUIRED	COLLECTION DATE
Syphilis Serology: VDRL FTA-ABS (FTA only for diagnostic problems)		Serum / CSF Serum	
Other Serology: Toxoplasmosis Rubella CMV Herpes Measles (Rubeola) Varicella Mumps Other, specify: _____		Serum Serum Serum Serum Serum Serum Serum Serum	
Agent Detection: Bordetella Pertussis Herpesvirus Agent Isolation: Streptococcus Group "A"		FA Slide FA Slide Throat Swab in Silica Gel	
Viral Isolation: (Specify Agent Suspected)		Throat Nasal Rectal Genital CSF Tissue Other _____	
Immunizations / Date		None Unknown MMR Varicella Influenza Adeno Other	
Contacts:		tick bite / date other / date	
Date of onset		Yes No	
Febrile			
Meningitis			
Respiratory			
Gastrointestinal			
Cardiovascular			
Paralytic			
Pregnant			
Other pertinent info:		(_____ weeks)	
Purpose of request:		diagnostic immune status antibody status other _____	

Date Received	Laboratory #	Technologist	Date Reported
---------------	--------------	--------------	---------------

Serum, 3 ml or Whole Blood, 6 ml	For antibody and/or immune status: submit a single serum For diagnostic determination: * by IgM tests - a single serum (a 2 nd serum may be requested later) by IgG tests -paired sera: <u>acute</u> phase collect within 7 days of onset of illness. <u>convalescent</u> phase collect 10 to 21 days later.	Ambient or refrigerated temperature
FA Slides for antigen detection	Slides/shippers provided by State Lab.	
Throat Swab (Group A Strep) in silica gel	Swab/mailler provided by State Lab	
Throat Washings	Use 5 - 10 ml sterile Hank's Balanced Salt Solution or sterile saline. 3 or 4 washings from the patient may be pooled in a sterile screw-cap jar. Seal tightly.	Specimens arriving within 24 hours of collection may be shipped refrigerated. If there will be longer storage or shipping times freezing is best. **
Throat Swabs, Rectal Swabs, Vaginal / cervical Swabs	A swab collection outfit is provided by the State Lab but any <u>Viral</u> transport medium is acceptable.	
Spinal Fluid	Submit in a sterile screw-cap tube and seal tightly.	
Feces	Place in a sterile container and seal tightly.	
Vesicle Fluid, or Pustule Crusts, Autopsy Tissues	Collect fluids on swabs and place in 1 ml of liquid (sterile Hank's or saline) in screw-cap container and seal tightly. Place each tissue in a separate sterile screw-cap container and seal tightly. Label each organ.	
Urine, fresh 10 ml	Place in a sterile screw-cap container and seal tightly.	Must be transported to the laboratory within 4 hours.

Specimen Submission: Select appropriate specimens for the clinical symptoms present, collect at proper intervals during illness, and handle as indicated below.

**Specimen Required
(without preservatives)**

Preparation

Shipping

* A significant rise in antibody ratios or titers may determine a recent infection. It is important that there is enough time between acute and convalescent sera to allow for antibody rise. A history of transfusion within 6 weeks of serum collection will invalidate serologic test results.

** When it is necessary to ship clinical samples frozen, use enough dry ice to last the trip. Seal the sample container with waterproof tape to protect the specimen from the CO₂. Wrap in absorbent material to minimize breakage and to contain any spills.

If the agent suspected is Respiratory Syncytial Virus (RSV) or Cytomegalovirus (CMV) the specimen should be refrigerated but do not freeze. It should arrive at the laboratory within 4 hours of collection.

A completed submission form #275 must be enclosed for each patient !
Dates of specimen collection, and when appropriate, date of onset of illness are required!

Rabies Agent Detection And/Or Isolation To Identify Infection In Animals

Persons Qualifying:

- All

Specimen:

- Animal head

Collection Kit (Rabies Kit) Furnished by State Lab Contains:

- Plastic Bucket with Lid-absorbent material to be provided by submitter
- Outer Cardboard Shipping Carton
- Manila envelope with:
 - Lab Form # 254- Rabies Examination
 - Address label # 254F
 - Instruction Sheet #254d
- Plastic Bags & tie fasteners

Cabinet for Health Services
Department for Public Health
Division of Laboratory Services
100 Sower Blvd Suite 204 Frankfort KY 40601

Lab 254A (Rev 8/99)

RABIES EXAMINATION

Address for hand delivery of specimens only!

Incident Information

Kind of animal: ☐ Dog ☐ Cat ☐ Fox ☐ Skunk ☐ Bat ☐ Raccoon ☐ Other _____

Was animal: ☐ Owned ☐ Stray

Was animal vaccinated? ☐ Yes, date ____/____/____ ☐ No
day year

Symptoms suggestive of Rabies? ☐ No ☐ Yes _____

County of incident: _____

Reason for Request

Person bitten? ☐ Yes, (name) _____ ☐ Animals exposed _____
(area of body) _____

Person: ☐ Scratched ☐ Licked ☐ Touched ☐ No known exposure

Other human exposure (specify) _____ ☐ Wildlife Survey _____
(county)

Specimen Information

Animal: ☐ Killed ☐ Died (Date: ____/____/____) Packed for shipment: ____/____/____
day mo year day mo year

Identification

Preference: (Must be a person's name) **Owner** if known, or **Person exposed**

Name: _____ Phone: ____/____/____

Address: _____

City: _____ ST: _____ ZIP: _____ County: _____

Submitting County Health Department: _____

City: _____ Phone: ____/____/____

If applicable: Vet Clinic, or Reference Lab (Ref # _____)

Name: _____ Address: _____

City: _____ ST: _____ ZIP: _____ Phone: ____/____/____

All below for DLS use ONLY

Id#:

Received:

Lab #:

###

#

Phone Record

Preliminary Report: _____
Date / time: _____

To: _____
By: _____

Confirmatory Report: _____
Date / time: _____

To: _____
By: _____

Lab 254B (Rev 8/99)

Human Immunodeficiency Virus (HIV) Serology

Patients Qualifying:

- Test available through AIDS counseling and testing sites or Sexually Transmitted Disease Clinics ONLY.

Specimen:

- Serum

Collection Kit (Specify Test) Furnished by State Lab Contains:

- Lab Form #197 Human Immunodeficiency Virus
- Mailing Label: #351 Serology
- Mailing container
- Red stoppered tube

Note: Numbered HIV Lab Form #197 and numbered stickers provided by STD Program.

<p style="text-align: center;"> Kentucky Public Health Laboratory 100 Sower Blvd., North Loading Dock, P.O. Box 2020 Frankfort, Kentucky 40602-2020 Phone: 502/564-4446 Fax: 502/564-7019 Thomas E. Maxson, Dr. P.H., Director </p> <p style="font-size: small;"> <i>Please complete a separate form for each specimen. Green copy may be retained by the submitter.</i> </p>	<h2 style="margin: 0;">Human Immunodeficiency Virus Serology</h2> <div style="display: flex; justify-content: space-around; align-items: center; margin-top: 20px;"> <div style="text-align: center;"> C. A <small>S</small> </div> <div style="border: 1px solid black; padding: 5px; text-align: center;"> A. K <small>D V</small> </div> </div>		
PATIENT INFORMATION:			
<div style="border-bottom: 1px solid black; margin-bottom: 5px;">Name (Last, First, MI)</div> <div style="display: flex; justify-content: space-between; border-bottom: 1px solid black; margin-bottom: 5px;"> Social Security # Sex Race Age Birthdate </div> <div style="border-bottom: 1px solid black; margin-bottom: 5px;">Home Address</div> <div style="display: flex; justify-content: space-between; border-bottom: 1px solid black; margin-bottom: 5px;"> City State Zip Code County </div>			
Send Report To:			
<div style="border-bottom: 1px solid black; margin-bottom: 5px;">Submitter</div> <div style="border-bottom: 1px solid black; margin-bottom: 5px;">Street Address (PO BOX)</div> <div style="display: flex; justify-content: space-between; border-bottom: 1px solid black;"> City State Zip Code </div>			
Specimen Information: Specimen type: <input type="checkbox"/> Serum <input type="checkbox"/> Whole Blood <input type="checkbox"/> Other _____ Date of Collection _____			
Program: Has patient been previously tested: <input type="checkbox"/> Yes <input type="checkbox"/> No If yes, when (date) _____: previous results: <input type="checkbox"/> Negative <input type="checkbox"/> Positive <input type="checkbox"/> Indeterminate			
Reason For Testing: (Mark One) <div style="display: flex; justify-content: space-between; margin-top: 5px;"> <div style="width: 45%;"> <input type="checkbox"/> Counseling-Testing Site Volunteer <input type="checkbox"/> Confidential <input type="checkbox"/> Anonymous <input type="checkbox"/> Adult & Child Health Clinic Patient <input type="checkbox"/> Symptoms suggest HIV Infection <input type="checkbox"/> Risk factors for HIV Infection </div> <div style="width: 45%;"> <input type="checkbox"/> TB Patient <input type="checkbox"/> STD Clinic <input type="checkbox"/> Person in Custody of Social Services <input type="checkbox"/> Needlestick Injury <input type="checkbox"/> Other (prior approval required) _____ </div> </div>			
<div style="text-align: center; border-top: 1px solid black; border-bottom: 1px solid black; margin: 5px 0;"> Laboratory Findings: </div> <div style="margin-bottom: 10px;"> Specimen Unsatisfactory: <div style="display: flex; justify-content: space-between;"> <div style="width: 30%;"> <input type="checkbox"/> Broken in transit <input type="checkbox"/> Insufficient quantity </div> <div style="width: 30%;"> <input type="checkbox"/> Chylous <input type="checkbox"/> Laboratory Accident </div> <div style="width: 30%;"> <input type="checkbox"/> Hemolyzed <input type="checkbox"/> Other _____ </div> </div> </div> <div> ELISA- Enzyme-Linked Immunosorbent Assay Test: <input type="checkbox"/> Non-reactive: No serological evidence of antibody to HIV-1 <input type="checkbox"/> Repeatedly reactive: Supplemental testing required </div> <div style="margin-top: 10px;"> Supplemental Test Performed: Western Blot <input type="checkbox"/> Non-reactive: No antibody to HIV-1 detected <input type="checkbox"/> Indeterminate: Testing inconclusive – Please submit an additional specimen in 6 weeks <input type="checkbox"/> Reactive: Antibody to HIV-1 detected </div>			
Date Received:	Laboratory Number:	Date Reported:	Technologist:

Please Use "L" Label or Fill in Completely

Prenatal Profile

Tests included in Prenatal Profile: ABO, Rh and antibody screen, Syphilis, HbsAg, Rubella

Patients Qualifying:

- Prenatal Patients

Specimen:

- Two full red-stoppered tubes of whole blood

Collection Kit (Prenatal Profile) Furnished by State Lab Contains:

- Two sterile red-stoppered tubes
- Lab Form: #212 Prenatal Profile
- Mailing Label: #359 Prenatal Profile
- Inner mailing container, absorbent material provided by submitter
- Outer mailing container

Notes or Comments: Use this profile only when ordering complete profile of tests. See lab form #213 for individual Syphilis, Rubella, Hepatitis, ABO, or Rh & antibodies testing.

Kentucky Public Health Laboratory
100 Sower Blvd., North Loading Dock,
P.O. Box 2020
Frankfort, Kentucky 40602-2020
Phone: 502/564-4446 Fax: 502/564-
7019

Prenatal Profile

Use this form for complete profile only,
see Lab Form 213 for individual tests.

(Please submit a completed Prenatal Profile Form
and two full 7 mL red stoppered tubes per patient.)

Patient Information (Please use L label or fill in completely):

Patient Name (Last, First, MI)

Patient I.D. # Sex Race Age DOB

Home Address

City State Zip County

Submitter Name Submitter Site Code

Weeks Pregnant Antepartum RhoGAM Date Date Collected

Prenatal Profile (ABO, Rh, and Antibodies, VDRL, HBsAg, Rubella) requires
two full 7 mL red-stoppered tubes.

FTA-ABS only performed when box below is checked.

☐ Do FTA-ABS if VDRL reactive. (Request only if patient has NO history of reactive FTA-ABS.)

Comments:

For Laboratory Use Only

White copy submitted with specimen • Yellow copy retained by submitter

HBsAg (Hepatitis B surface antigen), Anti-HBs (Antibody to HbsAg), or Anti-HBc (Antibody to HB core antigen)

Patients Qualifying:

- Prenatal patients, their contacts, and local health department employees (See Notes).

Specimen:

- Serum

Collection Kit (Hepatitis B) Furnished by State Lab Contains:

- Red-stoppered tube
- Lab form #213 Serodiagnosis
- Mailing label #351 Serology
- Inner mailing container, absorbent material provided by submitter
- Outer mailing container

Notes or Comments: Hepatitis B testing of local health department employees or patients other than prenatal patients should be approved by the Division of Epidemiology prior to testing.

ABO, Rh factor, Antibody Screen

Patients Qualifying:

- Prenatal Patients

Specimen:

- Whole Blood

Collection Kit (Prenatal Profile) Furnished by State Lab Contains:

- Red-stoppered tube
- Lab form #213 Serodiagnosis
- Mailing label #351 Serology
- Inner mailing container, absorbent material provided by submitter
- Outer mailing container

Syphilis Testing

Patients Qualifying:

- All local health department patients

Specimen:

- Serum

Collection Kit (Syphilis serology (single or double) Furnished by State Lab Contains:

- Red-stoppered tube
- Lab form #213 Serodiagnosis
- Mailing label #351 Serology
- Inner mailing container, absorbent material provided by submitter
- Outer mailing container

Special Serologies: *To order one of these tests,
mark "other serology" on reverse side and write in test name.*

Test:	Acceptable Specimen:
Candidiasis	Serum, 2 mL or Clotted Blood, 5-8 mL
Cryptococcosis	Serum, 2 mL or Clotted Blood, 5-8n mL or CSF, 1-2 mL
Hepatitis A IgM (Prior approval by the Division of Epidemiology required.)	Please call the Serology Section for instructions (502/ 564-4446).
Hepatitis C (Prior approval by the Division of Epidemiology required.)	Please call the Serology Section for instructions (502/ 564-4446).
Leptospirosis	Sera, 2 mL each; Acute and Convalescent
Rubella IgM (Prior approval by the Division of Epidemiology required.)	Please call the Serology Section for instructions (502/ 564-4446).
Sporotrichosis	Serum, 2 mL or Clotted Blood, 5-8 mL
Trichinosis	Serum, 2 mL or Clotted Blood, 5-8 mL

Note:

Hepatitis B testing of local health department patients other than prenatal patients and their contacts must be approved by the Division of Epidemiology prior to testing. Hepatitis B testing of local health department employees other than for determining immune status following immunization and in managing needlestick situations must also be approved by the Division of Epidemiology prior to testing.

<p style="text-align: center;"> Kentucky Public Health Laboratory 100 Sower Blvd., North Loading Dock, P.O. Box 2020 Frankfort, Kentucky 40602-2020 Phone: 502/564-4446 Fax: 502/564-7019 </p> <p style="font-size: small;"> <i>Please complete a separate form for each specimen. Yellow copy may be retained by the submitter.</i> </p>	<h2 style="margin: 0;">Serodiagnosis</h2> <p style="margin: 5px 0 0 0; font-size: small;">A double sided test order form</p>
PATIENT INFORMATION:	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> Name (Last, First, MI)	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> Social Security # Sex Race Age Birthdate </div>	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> Home Address	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> City State Zip Code County </div>	
Send Report To:	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> Submitter	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> Street Address (PO BOX)	
<div style="border-bottom: 1px solid black; height: 1.2em; margin-bottom: 5px;"></div> <div style="display: flex; justify-content: space-between; font-size: small;"> City State Zip Code </div>	
Specimen Information: Date of Collection _____ Specimen type: <input type="checkbox"/> Serum <input type="checkbox"/> Whole Blood <input type="checkbox"/> CSF <input type="checkbox"/> Other _____	
Purpose of Examination: <div style="display: flex; flex-wrap: wrap; font-size: small;"> <div style="width: 33%;"><input type="checkbox"/> Diagnostic</div> <div style="width: 33%;"><input type="checkbox"/> Pre-Hepatitis vaccine</div> <div style="width: 33%;"><input type="checkbox"/> Immune Status</div> <div style="width: 33%;"><input type="checkbox"/> Recheck Specimen</div> <div style="width: 33%;"><input type="checkbox"/> Post-Hepatitis vaccine</div> <div style="width: 33%;"><input type="checkbox"/> Prenatal _____ weeks pregnant</div> <div style="width: 33%;"><input type="checkbox"/> Treatment follow-up</div> <div style="width: 33%;"><input type="checkbox"/> Needlestick Injury</div> <div style="width: 33%;"><input type="checkbox"/> Other, specify _____</div> </div>	
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> Routine Examination Requested: <input type="checkbox"/> Rubella IgG <input type="checkbox"/> Syphilis testing </div> <div style="width: 45%;"> Hepatitis B (See note on reverse side) <input type="checkbox"/> HBsAg (Surface Antigen) <input type="checkbox"/> Anti-HBs (Antibody to HBsAg) <input type="checkbox"/> Anti-HBc (Antibody to HB Core Antigen) </div> </div> <div style="text-align: center; margin-top: 5px;"> Special Examinations <input type="checkbox"/> Other Serology, Specify _____ </div>	
Laboratory Findings:	

Please Use "L" Label or Fill in Completely

Dental Fluoride (Supplement Program)

Patients Qualifying:

- **Preschool children without a community fluoridated water supply.**

Specimen:

- **Sample of water supply**

Mailing label:

- **#505b Dental Fluoride**

Collection Kit Furnished by:

- **Dental Program 502/564-3246**

Notes and Comments: Instruction sheet included with collection kit.

Lab 505C
Rev. 1-99

Kentucky Public Health Laboratory
100 Sower Blvd., North Loading Dock,
P.O. Box 2020
Frankfort, Kentucky 40602-2020
Phone: 502/564-4446 Fax: 502/564-
7019

A.
B. Fluoride Test
For
Supplement
Program

(Please complete a separate form for each water supply.)

Name of Child(ren):

Sex:

DOB:

_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____

Home Address

City

State

Zip Code

Name of Parent or Guardian:

Send Report To:

Office/Clinic

Street Address (P.O. Box)

City

State

Zip Code

County

()
Phone Number

Specimen Information:

Water Supply: ☐ Well ☐ Cistern ☐ City ☐ Bottled Water
☐ Other, specify _____

Laboratory Findings:

____ . ____ (parts/million) µg/mL

Date Received:

Laboratory Number:

Date Reported:

Technologist:

Smear exam for Neisseria gonorrhoeae (GC)

Patients qualifying:

- All local health department patients

Specimen:

- Prepared direct smear from urethral exudate.

Collection Kit (Gonorrhea slide) Furnished by State Lab:

- Slide Mailer
- Lab Form # 219 Special Microbiology
- Mailing Label: Pre-addressed 3½ x 6½” envelope marked: GO, HAND CANCEL, FRAGILE HANDLE WITH CARE
- Instruction Sheet

Salmonella/Shigella species and other enteric pathogens

Patients Qualifying:

- All local health department patients.

Specimen:

- Stool specimen in “ENTERIC PRESERVATIVE” found in collection kit provided.

Collection Kit (Enteric pathogens) Furnished by State Lab Contains:

- Glass bottle, 1 oz., labeled “Enteric Pathogens”, containing a preservative
- Wood spoon for sampling
- Inner mailing container, with absorbent
- Outer mailing container
- Lab form #219 Special Microbiology
- Mailing label #353 Enteric Pathogens
- Instruction Sheet

Notes of Comments: If organisms other than Salmonella or Shigella sp. are suspected please call the Special Microbiology Lab at 502/564-4446 for further instruction or check the Reference Laboratory List.

Intestinal Parasites Examination (See Notes)

Patients Qualifying:

- All local health department patients

Specimen:

- Stool in 10% formalin or
- Double vial set of 10% formalin and PVA

Collection Kit (Intestinal parasites or Intestinal parasites w/PVA)

Furnished by State Lab Contains:

- Glass bottle, 1 oz., labeled
- “10% formalin, harmful if swallowed”
- Wooden spoon for sampling
- Inner mailing container, with absorbent
- Outer mailing container
- Lab form #219 Special Microbiology
- Mailing label #349a Parasitology
- Instruction sheet

Notes of Comments: All stool specimens submitted in 10% formalin are concentrated and examined for intestinal parasites. If the doctor requests the two-vial system, the second vial will contain PVA fixative. A trichrome stain and examination will be performed on stool specimens submitted in PVA fixative. The PVA fixative must be within the expiration date for valid testing. A larger outer mailing container will be necessary to accommodate the double vial system.

Pinworms

Patients Qualifying:

- All local health department patients.

Specimen:

- See instruction sheet provided; Perianal folds impression

Collection Kit (Pinworm Swube Tube) Furnished by State Lab Contains:

- Plastic adhesive paddle in a tube
- Outer styrofoam mailing container
- Lab form #219 Special Microbiology
- Mailing label #349a (Pre-addressed 3½ x 6½" envelope marked: PINWORM, HAND CANCEL, FRAGILE HANDLE WITH CARE)
- Instruction Sheet

<div>Kentucky Public Health Laboratory 100 Sower Blvd., North Loading Dock, P.O. Box 2020 Frankfort, Kentucky 40602-2020 Phone: 502/564-4446 Fax: 502/564-7019 Samuel B. Gregorio, Dr. P.H., Director</div> <div>Please complete a separate form for each specimen. Yellow copy may be retained by the submitter.</div>	<div>Special Microbiology</div>
PATIENT INFORMATION:	
Name (Last, First, MI)	
Social Security #	Sex Race Age DOB
Home Address	
City	State Zip Code County
Send Report To:	
Submitter	
Street Address (PO BOX)	
City	State Zip Code
Specimen Information:	
Purpose of Exam	<input type="checkbox"/> Clinical Specimen
Specimen Source	<input type="checkbox"/> Referred Culture
Date of Collection	Bloody Diarrhea <input type="checkbox"/> Ye <input type="checkbox"/>
Examination Requested: (Please mark one)	
<input type="checkbox"/> Smear Exam for GC	<input type="checkbox"/> Enteric Pathogens
<input type="checkbox"/> Direct Smear <input type="checkbox"/> Smear from Culture	<input type="checkbox"/> *Miscellaneous Bacterial Culture
<input type="checkbox"/> Culture Confirmation of Neisseria gonorrhoeae (GC)	<div>Organism Suspected:</div>
<input type="checkbox"/> Intestinal Parasites	
<input type="checkbox"/> Pinworm Prep	
<input type="checkbox"/> Other	
Other pertinent Medical Data: *Please complete this section when submitting Miscellaneous Bacterial Cultures	
FOR LABORATORY USE ONLY:	
Date Received:	Laboratory Number:

Please Use "L" Label or Fill in Completely

Test	Acceptable Specimen	Preservative/ Comments:
Smear Exam for <u>Neisseria gonorrhoeae</u> (GC)	Direct urethral exudate	None
Culture confirmation of <u>Neisseria gonorrhoeae</u>	Culture on applicable culture media	CO ₂ environment
Enteric Pathogens	1. Stool 2. Rectal Swabs 3. Culture on applicable culture media	Enteric Pathogens Kit Please call Special Bacteriology at 502/564-4446 for instructions.
Miscellaneous Bacterial Culture	Culture on applicable culture media	Please indicate any pertinent medical data, such as: clinical diagnosis; recent surgery/transplant; animal bites; diabetes, liver disease, etc.
Intestinal Parasites	Stool	10% Formalin
Miscellaneous Parasites	Specimens, such as arthropods, insects, or adult parasites	Please call Special Bacteriology at 502/564-4446 for instructions.
Pinworm Prep	See special collection instructions included in state kit.	Adhesive collection paddle, provided in state kit.

Mycobacteriology (TB) Smear & Culture

Patients Qualifying:

- All local health department patients

Specimen:

- Sputum

Collection Kit (TB sputum) Furnished by State Lab Contains:

- Sterile glass 1 oz. Vial with lid
- Inner mailing container, with absorbent
- Outer mailing container
- Lab form #207 MYCOBACTERIOLOGY SMEAR & CULTURE
- Mailing label #355 Pre-Paid, Mycob

Kentucky Health Services Laboratory 100 Sower Blvd., North Loading Dock, P.O. Box 2020 Frankfort, Kentucky 40602-2020 Phone: 502/564-4446 Fax: 502/564-7019 Samuel B. Gregorio, Dr. P.H., Director <i>Please complete a separate form for each specimen. Yellow copy may be retained by the submitter.</i>		Mycobacteriology Smear & Culture	
PATIENT INFORMATION:			
Name (Last, First, MI)			
Social Security #	Sex	Race	Age DOB
Home Address			
City	State	Zip Code	County
Send Report To:			
Submitter			
Street Address (PO BOX)			
City	State	Zip Code	
Requesting Physician (if other than submitter)			
Specimen Information:			
		Date of Collection _____	
<input type="checkbox"/> Clinical Specimen	<input type="checkbox"/> Referred Specimen		
<input type="checkbox"/> Sputum	Source: _____		
<input type="checkbox"/> Bronchial Washing			
<input type="checkbox"/> Gastric fluid	Hospital or Laboratory reference number		
<input type="checkbox"/> Urine	(if applicable _____)		
<input type="checkbox"/> CSF			
<input type="checkbox"/> Other, please specify _____			
Is the patient on anti-tuberculosis drugs? <input type="checkbox"/> Yes <input type="checkbox"/> No			
Laboratory Findings:			
Date Received:	Laboratory Number:	Date Reported:	Technologist:

Chlamydia trachomatis & Neisseria gonorrhoeae

Patients Qualifying:

- Family Planning, Prenatal, and STD Clinic Patients

Specimen:

- Cervical or Urethral Swab
- Specimen must be received in the lab within 7 days of collection. Transport at room temperature

Collection Kit (Chlamydia – GC (female) or Chlamydia – GC (male)) Furnished by State Lab Contains:

- Gen-Probe PACE 2 Specimen Collection Kit, specific for cervical or urethral site
- Inner mailing container
- Lab form #194 CHLAMYDIA TRACHOMATIS & NEISSERIA GONORRHOEAE by Nucleic Acid Probe
- Mailing label #194a Chlamydia and GC
- Outer mailer (single or multi size)

NOTE: Please discontinue use of the old single size styrofoam shippers previously used for the Chlamydia and use the double can system for 1-4 specimens or the multi-mailers for up to 17 specimens.

<p>Kentucky Public Health Laboratory 100 Sower Blvd., North Loading Dock, P.O. Box 2020 Frankfort, Kentucky 40602-2020 Phone: 502/564-4446 Fax: 502/564-7019</p>	<p>CHLAMYDIA TRACHOMATIS and NEISSERIA GONORRHOEAE by Nucleic Acid Probe</p>
<p><i>Yellow copy may be retained by the submitter.</i></p>	<p>D. See information on reverse side</p>
<p>PATIENT INFORMATION:</p>	
<p>Name (Last, First, MI) _____ (Codes defined on reverse side) 1 2 4 5 6 7</p>	
<p>Social Security # _____ Sex _____ Age _____ DOB _____ Race/Ethnicity (circle one)</p>	
<p>Home Address _____</p>	
<p>City _____ State _____ Zip Code _____ County _____</p>	
<p>Send Report To: _____</p>	
<p>_____ Health Department</p>	
<p>Street Address (PO BOX) _____</p>	
<p>City _____ State _____ Zip Code _____</p>	
<p>Reason For Testing: Did the patient present with Chlamydia/GC symptoms? <input type="checkbox"/> Yes <input type="checkbox"/> No Is the patient pregnant? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Unknown</p>	
<p>Mark one: <input type="checkbox"/> Volunteer/Medical Problem <input type="checkbox"/> Sex Partner Referral <input type="checkbox"/> Initial (Fam. Plan.) Visit <input type="checkbox"/> Other, please specify _____ <input type="checkbox"/> Revisit/Annual (Fam. Plan.) <input type="checkbox"/> Unknown/Undetermined _____ <input type="checkbox"/> Prenatal Visit <input type="checkbox"/> Cancer</p>	
<p>Specimen Information: Source (mark one): <input type="checkbox"/> Cervical <input type="checkbox"/> Vaginal, post-hysterectomy <input type="checkbox"/> Urethral <input type="checkbox"/> Other, specify _____</p>	
<p>Date of Collection _____ (dd-mmm-yy)</p>	
<p>~~~~~For Laboratory Use Only~~~~~</p>	
<p>Laboratory Findings:</p>	
<p><u>Chlamydia trachomatis</u> <input type="checkbox"/> Negative <input type="checkbox"/> Confirmed Positive <input type="checkbox"/> Presumptive Positive <input type="checkbox"/> Indeterminate <input type="checkbox"/> Unsatisfactory</p>	
<p><u>Neisseria gonorrhoeae</u> <input type="checkbox"/> Negative <input type="checkbox"/> Confirmed Positive <input type="checkbox"/> Presumptive Positive <input type="checkbox"/> Indeterminate <input type="checkbox"/> Unsatisfactory</p>	
<p>Remarks:</p>	
<p>Date and Time Received: _____ Laboratory Number: _____</p>	
<p>Date Reported: _____ Technologist: _____</p>	

Ethnicity: The Region IV Chlamydia project requires patient ethnicity data using the following codes:

- 1 White, Non Hispanic
- 2 Black, Non Hispanic
- 4 American Indian
- 5 Oriental
- 6 Hispanic White
- 7 Hispanic Black

Test: *Chlamydia trachomatis* & *Neisseria gonorrhoeae* by Nucleic Acid Probe

Specimen: Female Cervical or Male Urethral Swabs collected with the appropriate site specific Gen-Probe PACE 2 Specimen Collection Kit. For information concerning post-hysterectomy specimen collection, please call the Special Microbiology Section at 502/564-4446.

Specimen collection kits are supplied by this laboratory. Only use specimen collection kits within the stamped expiration date. Transport specimens at room temperature. Specimens need to be received in the laboratory within 7 days of collection.

Newborn Screening

Tests for the following are performed in the newborn screening laboratory:

- Phenylketonuria (PKU)
- Galactosemia
- Congenital Hypothyroidism (CH) [both T4 and TSH tests are performed]
- Sickle Cell

Patients Qualifying:

- All infants born in Kentucky.

Specimen:

- Whole capillary blood applied to the current lot number filter paper in the manner as described on the back of Lab Form #228.

Collection Kit (Newborn screening) Furnished by State Lab Contains:

- Lab Form #228 Newborn Screening.

MUST be current lot# of form and filter paper. Expiration Date indicates that the form is good through the month and year stated. For example, 04/03 indicates that the form is good through the last day of April 2003. Specimens collected after that date will be Rejected.

- Mailing Label: None Provided.

Cost:

- A charge of \$14.50 will be billed for those submitting an initial newborn screen. No charge will be billed for repeat specimens.

THE NEWBORN SCREENING FORMS, LAB FORM #228, REQUIRES VERTICAL STORAGE IN A COOL, DRY PLACE. DO NOT STORE IN PLASTIC BAGS.

Sickle cell Screening

Tests for the following are performed in the newborn screening laboratory:

- Sickle Cell

Patients Qualifying:

- Prenatal patients and older children.

Specimen:

- Whole capillary blood applied to the current lot number filter paper in the manner as described on the back of Lab Form #228.

Collection Kit (Newborn screening) Furnished by State Lab Contains:

Lab Form #228 Newborn Screening furnished by State laboratory.

MUST be current lot# of form and filter paper. Expiration Date indicates that the form is good through the month and year stated. For example, 04/03 indicates that the form is good through the last day of April 2003. Specimens collected after that date will be Rejected.

Mailing Label: None Provided.

THE NEWBORN SCREENING FORMS, LAB FORM #228, REQUIRES VERTICAL STORAGE IN A COOL, DRY PLACE. DO NOT STORE IN PLASTIC BAGS.

Form available from DLS

L-8844201 S&S® 903™ LOT W-001 Exp. 4/03

PHENYLKETONURIA (PKU), GALACTOSEMIA, CONGENITAL HYPOTHYROIDISM (CH), and SICKLE CELL SCREENING
 Samuel B. Gregorio, Dr. P.H., Director
 P.O. Box 2010 Frankfort, KY 40602
 Tel. # (502) 564-4448 Ext. 4434

Mother's First Name _____ Mother's Last Name _____
 Mother's SSN _____ County of Residence _____
 Street Address (PO Box) _____ State _____ Zip Code _____
 City _____
 Name of Infant's Physician _____ Infant's Physician's Phone Number _____
 Physician's Title _____ Street Address (PO Box) _____
 City _____ State _____ Zip Code _____
 Submitted by: _____
☐ Hospital ☐ Doctor ☐ Health Dept. ☐ Mobile
 Address _____ Ky _____
☐ First Specimen ☐ Repeat Specimen
☐ Newborn Screening (includes PKU, Galactosemia, T4/TSH and Sickle Cell)
☐ Sickle Cell Only (for older children and prenatal)

PRINT all information requested
 Infant's First Name _____ Last Name _____
 Date of Birth: Mo. _____ Day _____ Yr. _____ Time: AM _____ PM _____
☐ Male ☐ Female
 Ethnic Origin: ☐ White ☐ Black ☐ Other
☐ Hispanic ☐ Asian
☐ Premature ☐ Full Term
☐ Twins ☐ Triplets
☐ Antibiotics (ONLY IF GIVEN IN PAST 5 DAYS)
☐ Blood Transfusion (DATE _____)

NOTE: All infants tested before 48 hours of age must be retested before 3 weeks of age.
AGE OF INFANT WHEN SPECIMEN WAS DRAWN:
☐ 1-24 hours ☐ 25-47 hours ☐ 48-72 hours ☐ 73 hours-7 days ☐ 8 days-4 weeks ☐ Over 4 weeks
☐ Bottle ☐ Breast ☐ Both ☐ Feeding
 Birth Weight: _____ Grams
 Date of Birth: Mo. _____ Day _____ Yr. _____
 Date of Specimen: Mo. _____ Day _____ Yr. _____ Time: AM _____ PM _____
 Collector's Initial _____

200

S&S® 903™ LOT W-001 LAB 228 (REV. 2/2001)

Collection Instructions on Back of Form
 Submitter to retain pink copy.
 Fill circles completely using only one drop of blood per circle.

Front

NEWBORN SCREENING PROGRAM
PHENYLKETONURIA (PKU), GALACTOSEMIA, CH (CONGENITAL HYPOTHYROIDISM) and SICKLE CELL

- Obtain a specimen from each infant regardless of age, before the infant leaves the hospital. This excludes Neonates transferred to a higher level of care. Screen premature, ill or infants on parenteral feeding on the 7th day of life. Repeat screening is required for infants on antibiotics or parenteral feeding, or infants who received transfusions. Specific requirements for repeat screening are included in 902 KAR 4:030.
- It is recommended that specimens be collected prior to blood transfusion. The sickle cell test will be valid at this time.
- All infants tested before 48 hours of life MUST be retested prior to reaching 3 weeks of age for PKU, T4-TSH, and galactosemia.

INSTRUCTIONS FOR SPECIMEN COLLECTION

- DO NOT DETACH FILTER PAPER FROM FORM. DO NOT ALTER FORM.
- Cleanse the skin with an alcohol swab. Wipe off excess alcohol with dry sterile gauze.
- Puncture heel with sterile disposable lancet. Wipe away the first drop of blood with sterile gauze. Gently touch the filter paper against a large drop of blood. Blood spot should be large enough to soak through in ONE STEP. ALWAYS APPLY BLOOD TO ONE SIDE ONLY. NEVER APPLY ADDITIONAL BLOOD TO A FILLED CIRCLE.
- Allow blood specimen to AIR DRY THOROUGHLY on level non-absorbent open surface, such as a plastic-coated test tube rack for at least 3 hours. DO NOT HEAT, STACK, OR ALLOW BLOOD SPOTS TO TOUCH OTHER SURFACES DURING DRYING.
- SPECIMENS MUST BE MAILED WITHIN 24 HOURS OF COLLECTION.
- Check Sickle Cell only block when requesting a Sickle Cell test for an older infant or child.

- IT IS IMPERATIVE THAT ALL INFORMATION BE THOROUGHLY COMPLETED FOR ALL SPECIMENS SUBMITTED FOR TESTING.

S&S® 903™ LOT W-001 Exp. 4/03

S & S® 903™ LOT # W-001

Back

Cholesterol Screening (Must specify Total Cholesterol or Lipid Profile)

Tests included in Lipid Profile:

- Total Cholesterol, Triglyceride, HDL, LDL, VLDL

Patients Qualifying:

- Family Planning Patients / Chronic Disease Patients

Specimen:

- Serum, Lipid Profile requires fasting specimen.

Collection Kit (Cholesterol/Lipid screening) furnished by State Lab contains:

- Red stoppered tube
- Lab form #230 Clinical Chemistry
- Mailing Label #305 Clinical Chemistry
- Inner mailing container, absorbent material provided by submitter
- Outer mailing container

Notes/Comments: Please complete the risk Factors Box on the Submission Form for all Chronic Disease Patient testing.

Lithium Level

Patients Qualifying:

- Patient must have a Comprehensive Care physician's order prior to testing.

Specimen:

- Serum

Collection Kit (Lithium level) furnished by State Lab contains:

- Red stoppered tube
- Lab Form #230 Clinical Chemistry
- Mailing Label #305 Clinical Chemistry
- Inner Mailing Container, absorbent material provided by submitter.
- Outer Mailing Container

**Glucose Tolerance, Prenatal / Postpartum, One hour post 50
gram load Glucose**

Patients Qualifying:

- Prenatal / Postpartum patients

Specimen:

- Whole blood collected in sodium fluoride / potassium oxalate gray stoppered tube.

Collection Kit (Prenatal glucose tolerance) furnished by State Lab contains:

- Four gray stoppered tubes (one for each specimen)
- [Lab Form #230 Clinical Chemistry](#)
- [Mailing Label #305 Clinical Chemistry](#)
- Inner mailing container, absorbent material provided by submitter
- Outer Mailing Container or Small Styrofoam Box
- Cold pack

Notes / Comments: Please mail all timed tolerance specimens for one patient together with one completed submission form.

Fasting / Random Glucose

Patients Qualifying:

- Family Planning patients and Chronic Disease patients

Specimen:

- Whole blood collected in sodium fluoride / potassium oxalate gray stoppered tube.

Collection Kit (Prenatal glucose tolerance) furnished by State Lab contains:

- Four gray stoppered tubes (one for each specimen)
- Lab Form #230 Clinical Chemistry
- Mailing Label #305 Clinical Chemistry
- Inner mailing container, absorbent material provided by submitter
- Outer Mailing Container or Small Styrofoam Box
- Cold pack

TEST LIST	SPECIMEN REQUIRED
Total Cholesterol (only)	At least 2 mL of serum (4 mL of whole blood) required.
Lipid Profile Tests included: Total Cholesterol Triglyceride HDL LDL VLDL	At least 2 mL of serum (4 mL of whole blood) required. Specimen must be fasting for reliable results.
Prenatal/Postpartum Glucose Tolerance Test	At least 2 mL whole blood* per specimen is required. Complete one submission form per patient, mail all tubes together. Specimens must be mailed in an appropriate container with ice pack. *Sodium fluoride/potassium oxalate preservative is preferred, (gray topped tube).
Plasma Glucose	At least 2 mL whole blood* is required. Specimens must be mailed in an appropriate container with ice pack. *Sodium fluoride/ potassium oxalate preservative is preferred, (gray topped tube).

Clinical Chemistry

PATIENT INFORMATION:

Zip Code

Collection time

*Specimens must be mailed in an appropriate container with ice pack

☐ Sedentary Lifestyle